

3500927 Mop frame Clino® Synchro PRO Hybrid Mop frame for Cleanroom

Individualization
on request!

NEW



Product description

Autoclavable Mop frame for the cleanroom.

- Hygienic design
- Extremely robust locking mechanism
- Pretensioning of the sashes for optimum pressure distribution
- Simple mop holder
- Pretensioning of the sashes for optimum pressure distribution
- Contact-free mop discharge possible

**MADE IN
GERMANY**

ISO	3	4	5	6	7	8
GMP			A/B		C	D

Product specifications

Article name Article number	Mop frame Clino® Synchro PRO Hybrid 3500927
Dimensions Weight	40 x 10 cm 439 g
Opening symmetry	Symmetric
Wing material	Polyamide 6 Stainless steel 1.4404
Opening mechanism material	POM
Connection mechanism	Clino® Connect (Lock)
Connection Mechanism material	Polyamid 6 (Heat-stable special plastic)
Sterility	Can be autoclaved when open (121 °C for 20 minutes)
Care instructions	Please note the care instructions in the downloads.
Packing unit Selling unit	1 Pieces 1 PU
Quantity in carton	20 Pieces
Carton weight	8.78 kg
Cartons per pallet Pallet	no pallet shipment 120 x 80 x 220 cm
Packaging note	The size of the packaging, as well as the quantity in the box, may vary depending on the order quantity and the volume of other products ordered.
Customs tariff number	96039010
Country of origin	Germany
EAN-Number	4065237110645
Compatible accessories	Our sales department will be happy to advise you on suitable accessories

Compatible methods Cleaning area

- MopScoop®
- MopFloat®
- EasyMop®
- Flat wringer Ringo®



Industry Sectors

Hygiene

This classification is a non-binding recommendation. The user must check the individual suitability and compatibility for his individual area of application. The condition of the material must be checked regularly; damaged products must be replaced immediately. All information without guarantee. Subject to changes and deviations.

PFENNIG Reinigungstechnik GmbH | Heubachstraße 1 | 87471 Durach | Germany
www.pfennig-reinigungstechnik.com | info@pfennig-reinigungstechnik.com | +49 (0) 831 561 22 0

Status: 21.05.2025;

